

EC Sense EC Maintenance Option

12-lead Rest & Exercise ECG



CARDIOLEX


EASY TO USE ECG

**Streamline
your digital
workflow!**

EC Maintenance minimizes errors and downtime with step-by-step guides for preventive maintenance and simplified error logging – increases the efficiency for end users.



PROTOCOL OPTION




Preventive maintenance protocol




EASY TO USE

EC Maintenance is an option for EC Sense designed to streamline troubleshooting and preventive maintenance (PM). All Cardiolex products are characterized by its simple and intuitive use, which eliminates the risks for improper use. It is quick to get started with EC Maintenance, just select the chosen action under the menu items and start.




EASY PREVENTIVE MAINTENANCE WORKFLOW

The user can either perform troubleshooting or be helped through the various steps of the simple digital guide to perform preventive maintenance (PM) of the ECG system. The EC Maintenance control includes signal integrity tests, lead fail tests and function checks.



ERROR LOGGING

Troubleshooting logs with screenshots are easily generated by a keyboard shortcut from the device. The logs are saved to a file directory, USB stick or other media. In order to send log files and troubleshoot the system remotely, the ECG storage system EC Store and a basic EC Store support agreement are required.



PROTOCOL OPTION

The digital guide helps the user through the steps in the PM protocol and produces a final protocol in PDF. EC Maintenance has amongst other features a high sensitivity for detecting electrode errors and crosstalk between shield and leads, which is clearly shown in the report.



CARDIOLEX

EASY TO USE ECG

Cardiolex Medical AB develops and markets easy to use systems for ECG.

Our ECG management system and various ECG systems exist on 2/3 of Swedish hospitals, which means that the system has over 10 000 users in Sweden. We offer the complete solution within ECG, from the systems receiving the ECG signal from the patient to the ECG management systems, storing and managing generated data. Our products work for a wide range of users, from University hospitals to small clinics. Our focus is on creating ECG solutions that are user-friendly, reliable and secure.

Read more about our products on www.cardiolex.com.

Cardiolex Medical AB

Vretenvägen 6, SE-171 54 Solna

+46 (0) 8 564 820 50 · info@cardiolex.se

www.cardiolex.com